FLOCARE® BENCHMARK® NI TUBE

October 2020

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UK-IE-AU-NZ • BENCHMARK® NASO-INTESTINAL TUBE

PRODUCT DESCRIPTION

This is a medical device. The Flocare® Benchmark® Naso-Intestinal tube is a transnasal intestinal tube with a Benchmark® spiral tip design for the administration of enteral feeding.

The tube has the following characteristics:

Tube material: white full opaque polyurethane tubing

External diameter: Ch 8 or Ch 10

Length: 145 cm

Tip: bullet tip with 2 lateral outlets

Hydromer® coating on the tip

Benchmark® spiral (patented): 2,5 loops with a diameter of approx. 3cm and a length of approx. 23cm. These loops are straightened with a guide wire during insertion. After removal of the guide wire, the spiral eases transpyloric passage and allows optimal fixation in the small intestine.

Guide wire: A prelubricated metallic guide wire, partially introduced into the tube.

Connector: ENFit connector.

INDICATIONS

The Flocare® Benchmark® Naso-Intestinal tube is designed for transnasal feeding directly into the small intestine.

Indications for postpyloric tube feeding:

- Surgical patients
- Burn patients
- Patients with increased risk of aspiration, reflux and/or vomiting (often in early post-operative phase)
- Patients with limited gastric motility
- Patients with acute pancreatitis

Early postoperative intestinal tube feeding reduces the risk of atrophy of the small intestine and reduces the risk of bacterial translocation.

This tube is indicated for all patients who have a functioning intestinal tract with impaired stomach motility and/or increased risk of aspiration. This is frequently the case in the early postoperative period.

If there is stomach motility, the tube can be placed in the stomach (see instructions for use), the Benchmark® spiral will produce spontaneous pylorus passage within 8-12 hours.

If there is limited or no stomach motility, the tube can be passed through the pylorus with help of alternative techniques such as: endoscopy, or fluoroscopy. Alternatively the gastric motility can be improved by means of medication.

CONTRAINDICATIONS

For tube placement:

- Oesophageal injuries/haemorrhaging of the upper gastro-intestinal tract

For tube feeding:

- Intestinal absorption failure
- Paralytic or mechanical ileus
- Acute abdomen

INSTRUCTIONS FOR USE

Instructions for use for patients with (limited) gastric motility.

1. Advise the patient of the placement procedure. Fully insert the guide wire in the tube and ensure that it is firmly attached to the connector (fig.1).
2. Place the patient in a comfortable sitting or semi-recumbent position. To determine the length of Bengmark® tube to be inserted, measure the distance from the tip of the nose (A), to the ear (B), and then to the Xiphoid process of the sternum (C) (fig.2). Mark the tube at this point. Make 2 additional marks, at 25 cm and at 50 cm after the first mark.
3. Lubricate the tip with water, to facilitate insertion (fig.3).
4. Determine the preferred nostril and pass the tube along the bottom of the nostril (fig.4). When the tube reaches the throat, bend the patients' head forward and ask them to swallow as often as possible. Push the tube gently, never force it. Avoid inadvertent tube placement in the trachea! Continue until the first pre-measured mark is reached.
Warning: Do not connect the Bengmark® NI tube with the feeding set at this stage (fig 5). Correct placement must be confirmed before any fluid is passed through the naso-intestinal tube.
5. After checking the position of the tube (fig. 6, 7, 8), inject at least 20-50 ml H₂O (fig.9).
6. Pull ±25cm of guide wire out of the tube (fig.10). Then continue passing the tube until you reach the second mark (fig.11). Now pull the guide wire out of the tube completely.
7. Do not yet secure the tube to the nose, use the ear lobe as a temporary fixation point.
8. Within 8-12 hours gastric motility will have moved the tube through the pylorus (fig.12). Secure the tube once the third mark on the tube is at the patients' nose (fig. 13).
9. Before administrating any enteral nutrition, check if the tube is positioned correctly by X-ray (fig.14).

Giving something to drink or to eat (if possible) can stimulate the migration of the tube.

If required e.g. for paediatric use, the tube length can easily be reduced following tube insertion by cutting off the tube at the length required and inserting a new Flocare® feeding connector with the same Charrière size.

AFTER CARE

Flushing and Syringe use policy

- Use a syringe of at least 20ml for flushing the tube and always flush in a controlled manner.
- For medication delivery, always administer using an appropriate size syringe for the volume of medication but ensure delivery is in a controlled manner.
- Never use excessive force. If there is any resistance, stop the procedure and consult your clinician.

Check the position of the tube each time you change the tube feed, when in doubt whether the tube is in the correct position, and at least three times a day (fig. 6, 7, 8).

Flush the tube with 20-50 ml (fig. 15) of water before and after the administration of feed or medication and at least every 8 hours to prevent blockage of the tube.

It is recommended to place a new Flocare® Bengmark® Naso Intestinal tube every 6-8 weeks (fig.16). Replacement or removal of the product may differ from this recommendation, based on medical scrutiny, viability of the product and/or potential patient inconveniences.

It is recommended to always deliver feed in the small intestine using an enteral feeding pump to control the flow rate.

REMOVAL

The Benchmark® Naso Intestinal Tube can be removed by gently pulling the tube out of the patient's nose. Removal should always be performed under medical supervision.

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SINGLE PATIENT USE

UK-IE-AU-NZ • Do not reuse this product - Do not reuse, reprocess or re-sterilise this device. Reuse, reprocessing or re-sterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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