NUTRICIA
flocare[®]
GASTROSTOMY TUBE (G-TUBE)

UK-IE
GASTR

BE-NL
GASTR

DE-CH
GASTR

DK •
GASTR

ES •
SOND

FR-BE
SOND

FI •
GASTR

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IT-CH
SOND

NO •
GASTR

PT •
SOND

SE •
GASTR

EE •
GASTR

HU •
GASTR

LT •
GASTR

LV •
GASTR

PL •
ZGLE

TR •
GASTR

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
TR • KULLANIM KILAVUZU
GASTROSTOMİ TİPÜ

Manufacturer: Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp The Netherlands

CE
0344


STERILE EO


ENTERAL ONLY

200341_C

UK-IE

BE-NL

DE-CH

DK

ES

FR-BE

FI

GR

IT-CH

NO

PT

SE

EE

HU

LT

LV

PL

TR

www

- UK-IE-AU-NZ** Nutricia Ltd. TROWBRIDGE BA14 0XQ
Nutricia Ireland Ltd., Block 1, Deansgrange Business Park, Deansgrange, Co. DUBLIN
Nutricia Australia Pty Ltd, Bld D, Lvl 4, 12-24 Talavera Rd Macquarie Park NSW 2113
Nutricia Ltd 124 Favona Road, Favona, Auckland 2024, New Zealand
- BE-NL** N.V. Nutricia België,
Werkhuizenkaai 160, 1000 Brussel, België.
Nutricia Nederland B.V., Antwoordnummer 10230, 2700 VB Zoetermeer-NL,
Nutricia Medische Voedingsservice: 0800-022 33 22 (gratis).
- DE-CH-AT-LU** Nutricia GmbH, D-91052 ERLANGEN, Deutschland
Nutricia GmbH, A-1120 WIEN, Österreich
Nutricia SA, CH-8005 ZÜRICH, Schweiz - Produkt- und Ernährungsberatung: 00800 700 500 00
- DK** Nutricia A/S, Rørmosevej 2A, 3450 ALLERØD
- ES** Nutricia SRL, C/ Torrelaguna 77, 6ª plta. CP 28043, Madrid, España. www.nutricia.es
Tel. atención al cliente: (+34) 900 211 088.
- FR-BE-CH-LU** NUTRICIA Nutrition Clinique, 150 Bd Victor Hugo,
93406 Saint-Ouen Cedex, France. www.nutricia.fr
N.V. Nutricia België, 160 Quai des Usines, 1000 Bruxelles, Belgique.
Pour plus d'information : 070/222.307.
Nutricia S.A. Suisse, CH-8005 ZÜRICH.
- FI** Nutricia Medical Oy, Linnankatu 26A, 20100 TURKU
- GR** Numil Hellas AE 17 χλμ. Αθηνών Λαμίας & Καλαμάτας 2, 14564 Κηφισιά
Γραμμή φροντίδας:800 11 68600,210 6248 547
- IT-CH** Nutricia Italia SpA, via C.Farini 41 - 20159 Milano - Numero Verde: 800822096
Nutricia SA Svizzera, CH-8005 ZÜRICH
- NO** Nutricia Norge AS, Drammensveien 123. 0277 Oslo
- PT** Distribuído em Portugal por: Nutricia Advanced Medical Nutrition – Unipessoal, Lda.
Avenida D. João II nº41, 3º piso Torre Zen, 1990-084 LISBOA
- SE** Nutricia Nordica AB, Box 3142, 169 03 Solna
- EE** NutriMedical OÜ, TALLINN 11313, Eesti
- HU** Forgalmazza: Numil Kft. 1134. BUDAPEST, Váci út 35. www.nutricia.hu
- LT** Nutricia SIA, J. Gustava Zemgala gatve 74, Ryga, LV-1039, Latvija tel. +37052300377
- LV** Nutricia SIA, Gustava Zemgala gatve 74, Riga, LV-1039, Latvija
- PL** Nutricia Polska Sp. z o.o., ul. Bobrowiecka 8, 00-728 WARSZAWA
- TR** Numil Gıda Ürünleri San. ve Tic. A.Ş. Yapı Kredi Plaza B Blok Kat:8-9
Cömert Sok. Levent Mah. Beşiktaş- İSTANBUL

www.nutriciaflocare.com

N
3

00

s

:77

UK-IE-AU-NZ • GASTROSTOMY TUBE

PRODUCT DESCRIPTION

The Flocare® Gastrostomy tube (G-tube) is a silicone balloon catheter used as a replacement of an existing gastrostomy tube (PEG, G-tube or button). The tube is made of transparent silicone and contains:

- a balloon inflation port indicating the advised balloon filling volume (A)
- a feeding port (ENFit connection with closure cap) (B)
- cm markings on the tube
- a silicone external fixation disc (C)
- a silicone internal retention balloon; advised filling volume is dependent on Charrière size of the tube and indicated on the balloon inflation port (A)
- a quick release clamp; to prevent backflow of nutrition or other gastric content (D)
- a radiopaque mark near the entrance of the balloon enabling x-ray visibility (E)

The Flocare® G-tube is available in 5 Charrière sizes (10, 14, 16, 18 and 20)

CONTENTS

Each customer packaging box, contains two Flocare® Gastrostomy tubes, which are individually packed in a blister.

INDICATIONS

- Long term tube feeding exceeding the life time of the initially placed gastrostomy tube.
- Stomach decompression.

CONTRA-INDICATIONS

It is contra-indicated to place or replace a G-tube in case of:

- Non-established or infected stoma site,
- Ascites,
- Peritoneal carcinoma.

INSTRUCTIONS FOR USE

1. Position the patient onto their back (fig.1). Check the Charrière size and length of the existing tube. Apply anaesthetising gel around the tube. Place a gauze under the external fixation disc to prevent leakage of gastric liquid directly onto the skin.
2. Take the new Flocare® G-tube of the same Charrière size as the tube in situ. Check the integrity and state of the balloon by injecting the indicated amount of sterile water in the balloon inflation port (A) (fig. 2 and 3).
3. Roll the balloon, while wearing sterile gloves, between the thumb and the index finger to ensure it is centred on the tube. Deflate the balloon (fig.4).
4. Close the cap of the feeding port (B) (fig.5).

Remove the tube as per relevant instructions. Clean the exit of the stoma with dry gauze. Apply a water-soluble lubricant or an anaesthetising gel on the tip of the new tube, keeping the outlet of the tube tip clear.

5. Introduce the new tube as far as possible into the stoma (fig.6).
6. Keep the tube in place and confirm the position of the tip by aspirating gastric contents via the feeding port (B) and by measuring the gastric pH. The tip of the tube is situated in the stomach when the pH is lower than 5.5 (fig.7). When there is any doubt, confirm the position by means of an endoscope or X-ray. When opening the feeding port, the quick release clamp (D) prevents backflow of nutrition or other gastric content. It is recommended to open the quick release clamp (D) when the feeding port is closed.
7. Inflate the balloon via the balloon inflation port (A), with the quantity of sterile water indicated on the tube (see fig. 8). Never exceed the volume capacity of the balloon. Pull the tube slightly upwards until the length of the tube coming out of the stoma is identical to the length of the previous G-tube in situ (if relevant) (fig.9). Note the Charrière, length of tube introduced and the amount of water injected in the balloon, in the nursing and medical records.
8. Flush the tube with 20-50 ml of water via the feeding port of the tube (B) to avoid an obstruction of the tube (fig.10).
9. Slide the external fixation disc (C) approximately 5 mm from the skin and fix the tube in the 90°-angled disc (fig.11).

MAINTENANCE

Inspect the stoma-site daily for any redness and swelling. Disinfect the skin, when necessary. As soon as the stoma is completely formed it is sufficient to wash, rinse and dry the skin thoroughly.

Daily Care: remove the tube out of the external retention disc and gently introduce the tube further into the stoma. Turn the tube 180° on its axis, pull the tube upwards to its original position and fix the tube back in the external fixation disc (fig.12).

Check the position of the tube by checking the position of the external fixation disc or, if possible, by measuring the pH value of aspirated gastric content. pH value below 5.5 confirms a correct position of the tube in the stomach. Do this each time you connect a new packaging of tube nutrition or if there is any doubt whether the tube is in the correct position (fig 13).

Flush the G-tube with 20-50 ml of water before and after the administration of feed or medication and at least every 8 hours to prevent tube blockage (fig.14).

Check on a weekly basis the volume of sterile water in the balloon. If necessary fill up or replace the volume to the same amount as documented in the medical or nursing records.

Never administer tube feeding or medication into the balloon inflation port (A).

DURATION OF USE AND REMOVAL

The Flocare® G-tube is designed to have a life span of several months. The life span of the tube varies according several factors which include medication, gastric pH, patient motility, duration of use and tube care.

Replacement and removal of the G-tube depends upon the viability of the product and/or should be performed under medical instruction.

To remove the G-tube, simply deflate the balloon with a syringe and gently pull out the tube until the tip is out of the stoma; force should not be used.

If the balloon cannot be deflated, clean the deflation port with warm water and retry. If this is unsuccessful, it is recommended that any further attempt to remove the tube be carried out under medical instruction, preferably in a hospital setting.

It is not recommended to:

1. Fill the balloon to bursting as this can cause trauma
2. Cut the G-tube at skin level as this can cause obstruction
3. Empty the balloon via the tube with a needle or syringe as this cannot be performed with certainty and poses significant risk to the patient

In case of accidental removal: immediately introduce a new G-tube (of the same Charrière) to keep the fistula open. Inflation of the balloon may only be done under medical supervision or endoscopic control to prevent damage of the fistula.

Flushing and Syringe use policy

- Use a syringe of at least 20ml for flushing the tube and always flush in a controlled manner.
- For medication delivery, always administer using an appropriate size syringe for the volume of medication but ensure delivery is in a controlled manner.
- Never use excessive force. If there is any resistance, stop the procedure and consult your clinician.

- It is recommended to open the quick release clamp (D) when the feeding port is closed as a continuously closed clamp can damage the tube and material integrity.
- Material defects and tube occlusion can be largely ruled out if the Flocare® G-tube is adequately flushed before and after every administration of nutrition and/or medication
- Use liquid medication where possible
- Avoid using acidic substances such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging

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Gen
ved
re-ste
kryds
en p
Efter
eller

SINGLE PATIENT USE

UK-IE-AU-NZ • Do not reuse this product - Do not reuse, reprocess or re-sterilise this device. Reuse, reprocessing or re-sterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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