NUTRICIA
flocare®

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EE • BENCHMARK® PEG/J for use with FLOCARE® PEG CH 18

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UK-IE-AU-NZ • BENCHMARK® PEG/J for use with FLOCARE® PEG CH 18

This is a medical device. Contains loose parts which may be dangerous to young children / DEHP-free / Latex-free

Manufacturer: Nutricia Medical Devices BV, Taurusavenue 167, 2132 LS Hoofddorp, The Netherlands

Production location: Nutricia Pharmaceutical Wuxi Co., Ltd., No. 17 Xinming Road, Wuxi High-tech Development Zone, P.R. China

STERILE EO

ENTERAL ONLY

NUTRICIA
flocare[®]

FLOCARE[®] BENCHMARK[®] PEG/J for use with FLOCARE[®] PEG CH 18

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UK-IE-AU-NZ • BENCHMARK® PEG/J for use with FLO CARE® PEG CH 18

PRODUCT DESCRIPTION

The Flo care® Benchmark® PEG/J is a jejunal feeding tube which should be used with a Flo care® PEG CH 18. Nutrition/medication can be administered using an ENFit feeding set.

Tube material:	blue fully opaque polyurethane tubing (A)
External diameter:	CH 9
Length:	105 cm
Tip:	antenna-tip with 2 lateral outlets (B). Benchmark®-spiral (patented): 2,5 loops with a diameter of approx. 3 cm and a length of approx. 23 cm. These loops are straightened with a guide wire during insertion. After removal of the guide wire, the spiral allows optimal stable positioning in the small intestine. A black square shape marking printed on the tube about 25 cm from the tip indicates to which point the spiral will roll up when the guide wire is removed.
Guide wire:	a prelubricated metallic guide wire, partially inserted into the tube (C).
ENFit Connector:	consisting of lilac (D) and turquoise (E) parts, to enable connection of the jejunal tube with the Y-connector.
Y-Connector (red):	to allow simultaneous aspiration of gastric contents and feeding in the small intestine. The gastric leg (F) contains a funnel connector allowing connection of a step connector with a drainage bag. The jejunal leg (G) is pre-assembled on the jejunal tube. The Y-connector contains an inner ring (H) to attach the connector to the PEG, as well as a separate red ring (I) to fixate the connection. The turquoise ring (J) on top of the Y-connector enables connection with the jejunal tube.
Spare Connector:	to ensure optimal connection after tube placement; consisting of white, transparent (K) and turquoise (L) components which can be disconnected separately.
Distance marks:	every 10 centimeters.

INDICATIONS

The Flo care® Benchmark® PEG/J CH 9 in combination with the Flo care® PEG CH 18 is suitable for those patients requiring enteral tube feeding directly in the jejunum or duodenum via a gastrostomy with a jejunal feeding tube, potentially together with simultaneous gastric drainage.

Indications for postpyloric tube feeding:

- Surgical patients
- Burn patients
- Patients with increased risk of aspiration or reflux and vomiting (early postoperative phase)
- Patients with limited gastric motility
- Patients with acute pancreatitis

This group of patients typically have a functioning intestinal tract with impaired stomach motility and/or increased risk of aspiration, e.g. early postoperative feeding. Early postoperative intestinal tube feeding reduces the risk of atrophy of the small intestine and reduces the risk of bacterial translocation.

CONTRAINDICATIONS

- For PEG placement:
- Severe coagulation disorders (e.g. haemophilia)
 - Severe ascites and severe gastritis
 - Ileus
 - Anorexia nervosa and severe psychoses
 - Acute pancreatitis and acute peritonitis
 - When it is impossible to bring the anterior part of the stomach against the abdominal wall (subtotal gastrectomy, ascites, hepatomegaly, etc.)
 - Marked peritoneal carcinomatosis
 - Extensive tumor infiltration in the area of the puncture site
- For postpyloric feeding:
- Intestinal absorption failure
 - Paralytic ileus
 - Acute abdomen

More detailed information can be found in the Instructions For Use for the Flocare® PEG.

INSTRUCTIONS FOR USE

Preparation of the Flocare® PEG CH 18

After initial placement of the Flocare® PEG CH 18, with the endoscope still in stomach:

- Cut-off PEG tube at 35-cm marking.
- Pull external retention disc over the tube, but do not fix the tube in the disc yet.
- Place quick release clamp on the tube.
- Lubricate the inner-lumen of the PEG tube with approx. 10 ml water, using a syringe.

A. Endoscopic placement:

1. Insert the guide wire completely into the jejunal tube. Withdraw the guide wire about 10-20 cm, and then push the guide wire until it reaches the end of the tube. Attach the handle of the guide wire (C) firmly to the ENFit Connector (D).
2. Lubricate the stretched jejunal tube with water to facilitate intubation.
3. Push the antenna tip (B) of the jejunal tube through the PEG CH 18 until the tip is visible in the stomach. If necessary add some extra water through the PEG CH 18 tube.
4. Grasp the antenna-tip with the endoscope forceps (fig. 1) and with help of the endoscope guide the tube through the pylorus. The tube tip should be positioned as far down the small intestine as possible, preferably beyond the Treitz ligament (fig. 2).
5. Keep the tube tip in place with the forceps, whilst withdrawing the endoscope (withdrawal of the endoscope approximately 30 cm is possible) (fig. 3).
6. Whilst the tube is still held by the forceps, withdraw the guide wire about 25 cm from the tube. Now the spiral will begin to form. With the forceps still in place, move the blue jejunal tube forward further into the small intestine (fig. 4).
7. Finally remove the endoscope, the forceps and the guide wire cautiously (fig. 5) being careful not to move the tube backwards at the same time.
8. Detach the ENFit connector from the Y-connector by unscrewing the turquoise ring (J) on top of the Y-connector (fig. 6).
9. Unscrew the separate red ring (I) from the Y-connector and move it over the PEG CH18 (fig. 7).
10. Connect the Y-connector to the PEG CH18 by fitting the red inner shaft (H) tightly into the PEG tube. Screw the separate red ring (I) into the Y-connector to fixate the Y-connector on the PEG tube (fig. 8)
11. Cut the ENFit connector from the jejunal tube (5 cm above the Y-connector) (fig. 9). If necessary, check the position of the tip of the tube endoscopically, to rule out any dislocation at this time.

12. Unscrew the turquoise ring (L) of the spare connector from the transparent part (K). Slide the turquoise ring over the blue jejunal tube. Fix the jejunal tube over the metallic pin inside the transparent part of the spare connector. Make sure the tube is completely moved over the metallic pin. Connect the turquoise and transparent parts of the spare connector by screwing in the turquoise part. Afterwards, connect the complete assembled connector to the Y-connector, by screwing the turquoise ring of the Y-connector (J). (fig. 10 and fig. 11).
13. Fix the PEG CH 18 tubing in the external fixation disc. Ensure the correct positioning of the external fixation disc is maintained by placing the blue safety clamp immediately behind the disc.
14. Check that the jejunal tube works by flushing with approx. 20 ml water.
15. Check if the jejunal tube is in the correct position with an X-ray. Note the length of the remaining part of the intestinal tube.

B. Fluoroscopic placement:

Follow point 1-3 of the instructions mentioned under A.

4. Check that the tube tip is in the stomach with fluoroscopy and direct the tip towards the pylorus. Withdraw the guide wire approximately 2-3 cm which will cause the tube to bend. This can facilitate finding the pylorus and help inserting the tube-tip through it. The safety tip design of the Bengmark® PEG/J tube enables manipulation with the guide wire in situ without the risk of lateral guide wire penetration, which could cause tissue damage. After pylorus penetration the guide wire should be brought back in the original position.
5. Insert the Bengmark® PEG/J tube with the guide wire as far as possible in the small intestine, preferably beyond the Treitz ligament.
6. Once the tube is in place, withdraw the guide wire.
7. Follow point 8-15 of the instructions mentioned under A.

AFTER CARE

Verify status of the tube by checking the position and state of the external fixation disc and the tube itself, each time you change the tube feed, if in doubt whether the tube is in the right position, and at least three times a day.

Flush the tube before and after feed and/or medication administration, and at least every 8 hours with 20-50 ml of water to prevent tube blockage (see figure 13). The Flocare® Bengmark® PEG/J is used solely in combination with the Flocare® PEG CH 18, which in turn implies similar maintenance practices. Ensure you follow the maintenance instructions for the PEG carefully, also when used in combination with the Bengmark® PEG/J.

Be aware that maintenance with the PEG/J-PEG combination differs in one area (see below). Within the PEG instructions for use the following is mentioned: Release the tube out of the external fixation disc, turn the tube daily 180° on its axis, and move tube up and down in the fully healed stoma (min 15 mm). When the PEG is combined with a Bengmark® PEG/J, the tube should NOT be turned on its axis, but only moved up and down in the fully healed stoma (min 15 mm). Turning the tube is NOT carried out so as to avoid twisting of the Bengmark® PEG/J CH 9 tube, leading to possible dislocation of the PEG/J Bengmark® coil. Always return the tube to its original position with the PEG's blue safety clamp behind the external retention disk. (see figure 12).

Disinfectants such as Povidone-Iodine (PVP-I; e.g.: Iso-Betadine®, Braunol®) and Octenidindihydrochlorid-Phenoxyethanol (e.g.: Octenisept®) should not be used because repeated exposure to this disinfectant might negatively effect the physical/mechanical properties of the tube. It is recommended to use disinfectants based on Polyhexanid (e.g. Prontosan®)

**Flushing and Syringe use policy**

- Use a syringe of at least 20ml for flushing the tube and always flush in a controlled manner.
- For medication delivery, always administer using an appropriate size syringe for the volume of medication but ensure delivery is in a controlled manner.
- Never use excessive force. If there is any resistance, stop the procedure and consult your clinician.

REMOVAL

It is recommended to place a new Flocare® Benchmark® PEG/J tube every 6 weeks. Replacement or removal of the product may differ from this recommendation, based on medical scrutiny, viability of the product and/or potential patient inconveniences.

It is recommended to deliver tube feed using an enteral feeding pump, to control the flow rate.

SINGLE PATIENT USE

UK-IE-AU-NZ • Do not reuse this product - Do not reuse, reprocess or re-sterilise this device. Reuse, reprocessing or re-sterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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